



MALES

STANDARD URINE SPECIMEN COLLECTION

1. Wash hands with soap.
2. Open the container **WITHOUT TOUCHING THE INSIDE OF THE CUP**.
3. Use the towelette to clean the end of the penis. Discard and clean again with the second towelette. If not circumcised, hold foreskin back before cleansing.
4. Begin urinating into the toilet and while continuing to urinate, bring the cup into the stream and continue to urinate until cup is at least $\frac{1}{2}$ full.
5. Finish voiding into the toilet.
6. Replace lid on cup. Tighten securely.
7. Give the specimen to a technician or place in the designated area.

PACIENTES VARONES

RECOLLECTION DE ORINA

1. Lavesé las manos con jabón.
2. Quite la tapa del envase sin tocar el interior.
3. Limpiese la punta del pene con la toallita. Repitalo con una toallita más.
4. Empiece a orinar en el baso del baño y sin dejar de orinar colecte la orina en el envase llenandolo hasta la mitad.
5. Termine de orinar en el baso del baño.
6. Tape el envase de nuevo.
7. Dé le la muestra a un tecnico o pongalo en el lugar designado.