



Chronic Kidney Disease Test Menu

Name (s)	Test Codes	Specimen Requirements	Clinical Application
Estimated Glomerular Filtration Rate (eGFR)	N/A	eGFR is calculated from the Creatinine result anytime a Creatinine is ordered, whether individually or within a panel or profile. Please provide patient age and gender.	Following KDOQI guidelines, two consecutive eGFR calculations within a three month period will assist physicians in staging of Chronic Kidney Disease.
Phosphorus	1012	2 mL REFRIGERATED serum (1 mL minimum) from a serum separator tube or a plain red-top tube	Useful in monitoring treatment of patients with Chronic Kidney Disease.
PTH Intact & Calcium	102846	1 mL FROZEN serum (0.5 mL minimum) from a serum separator tube AND 2 mL REFRIGERATED serum (1 mL minimum) from a serum separator tube or a plain red-top tube	Useful in determining secondary hyperparathyroidism.
Vitamin D 1, 25-Dihydroxy	8115	3 mL FROZEN serum (1.1 mL minimum) from a plain red-top tube or serum separator tube	Useful in determining treatment of patients with Chronic Kidney Disease.
Vitamin D, 25-Hydroxy – LC/MS/MS	902068	0.25 mL ROOM TEMPERATURE serum (0.15 mL minimum) from a plain red-top tube or serum separator tube	Useful in determining treatment of patients with Chronic Kidney Disease.

For results or questions contact our Client Services Department:
Phoenix Location at 602.685.5050 or 800.766.6721, extension 45050 – toll free;
Tucson Location at 520.886.8101 or 800.266.8101 – toll free.
For further information on any of the testing on this menu, please call the Sales & Marketing Department:
602.685.5285 or 800.766.6721, extension 45285 – toll free.

